MADE IN ITALY

**MAGIC CARE®**

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## Mod. P0406EM F400

### Aerosol therapy unit

Congratulations on your purchase and thank you for choosing our unit. Our goal is to fully satisfy consumers by offering them cutting-edge systems for the treatment of respiratory tract ailments.

Carefully read these instructions and keep them in a safe place for future reference.

Use the unit only as described in this instruction manual. **This is a medical device for home use designed to administer nebulizer medications prescribed or recommended by your physician.**

Please visit [www.flaemnuova.it](http://www.flaemnuova.it) to see the entire range of Flaem products.

### THE STANDARD COMPONENTS AND ACCESSORIES OF THE UNIT INCLUDE:

#### **A) - Aerosoltherapy unit (main unit)**

- A1 - On/off switch
- A2 - Air outlet port
- A3 - Air filter
- A4 - Nebuliser holder
- A5 - Socket for power cord
- A6 - Power cord

#### **B) - Connection tube (main unit / nebuliser)**

#### **C) - Accessories**

- C1 - RF6 basic<sup>2</sup> nebuliser
  - C1.1 - Lower part
  - C1.2 - Complete nozzle
  - C1.3 - Upper part
  - C1.4 - Cap
- C2 - Mouthpiece
- C3 - Adult mask (PVC)
- C4 - Child mask (PVC)

#### **D) - Carrying bag**

### **⚠ IMPORTANT SAFEGUARDS**

- Before using the unit for the first time, and periodically throughout its life, check the power supply cord for any damage; if you detect any damage, do not plug it in. Take the appliance immediately to an authorised service centre or your dealer.
- The average lifespan of accessories is 1 year; however, we recommend you replace the nebulizer every 6 months in intensive use (or earlier, if the nebulizer is blocked) to always ensure maximum therapeutic effect.
- Do not allow to be used as a toy. Close attention is necessary when used by or near children.
- Children and disabled persons should always use the unit under strict supervision of an adult who has read this manual.
- Several parts of the device are small enough to be swallowed by children; therefore keep the device out of reach of children.
- Do not use the tubes and the cords of the device for purposes other than those specified as they may cause strangulation hazard; utmost attention must be paid to children and persons with special needs; usually, such persons are not capable of assessing the hazard correctly.
- The unit should not be used in the presence of anaesthetic mixture inflammable with air, oxygen or nitrous oxide.
- Keep the cord away from heated surfaces.
- Keep the power supply cord away from animals (for example rodents) as they may damage the insulation of the power supply cord.
- Do not handle the unit with wet hands. Do not use the device in humid environments (e.g. while taking a bath or shower). Do not immerse the device in water; should this accidentally occur, unplug the device immediately. Do not remove or touch the device immersed in water before pulling out the plug. Take the device immediately to an authorized service center or to your local dealer.
- Only use the device in dust-free environments; otherwise, therapy may be compromised.
- The casing of this unit is not protected against the penetration of liquids. Do not wash the device under running water or submerge in water. Keep away from splashing water or other liquids.
- Do not expose the device to high temperatures.
- Do not put the device near heat sources, direct sunlight, or in hot environments.
- Do not obstruct or insert any objects in the filter or in its housing.
- Do not obstruct the air vents placed on both sides of the unit.
- During use, always place the unit on a firm surface clear of obstacles.
- Check that there is no material obstructing the air vents before each use.
- Do not place any objects inside the air vents.
- The average duration provided for compressor ranges is: F400: 400 hours, F700: 700 hours, F1000: 1000

hours, F2000: 2000 hours.

- Repairs must be performed by authorized personnel only. Any unauthorized repairs will void the warranty and may pose a safety hazard for the user.


**WARNING:** Do not modify this equipment without authorization of the manufacturer

- The Manufacturer, the Retailer and the Importer are responsible for the safety, reliability and performance of the device only if: a) the device is used in accordance with the operating instructions b) the electrical system of the premises where the device is used is in compliance with current laws.
- Interactions: The materials used that come into contact with medications have been tested with a wide range of medications. However, due to the variety and continuous evolution of medication, the possibility of a chemical interaction cannot be excluded. We recommend that once the medication has been opened you should use it as soon as possible and avoid prolonged contact of the medication with the nebulizer.
- The manufacturer should be contacted for reporting problems and/or unexpected events related to device operation.


## INSTRUCTIONS FOR USE

**Before each use, the nebuliser and accessories should be cleaned according to the instructions in the section "CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION". A personal use of the nebuliser and its accessories is recommended in order to avoid any risks of infection.**

1. Plug the power cord (A6) into the socket (A5) at the rear of the unit and plug the power cord into a power socket corresponding to the voltage of the unit. This must be positioned so that it is not difficult to make the disconnection from the mains.
2. Open the nebuliser by turning the upper part (C1.3) counterclockwise.
3. Pour the medication prescribed by the physician into the lower part (C1.1). Close the nebuliser by turning the upper part (C1.3) clockwise.
4. Connect the accessories as shown in the "Assembly diagram". The cap (C1.4) provided with the nebuliser RF6 basic<sup>2</sup> makes it possible to deliver the medication in 2 modes:

**"standard mode"**; for optimum delivery and minimum waste of medication

**"high speed mode"**; for a quicker treatment

If you use the face mask, fit it to the face as shown in the picture

5. Start the unit by means of switch (A1) and breathe in and out deeply.
6. When the treatment is finished, turn off and unplug the unit.
7. If after the therapy session moisture accumulates in the tube (B), remove the tube from the nebulizer and dry it using the compressor air flow; this operation will prevent mold from forming inside the tube.

## CLEANING, SANITIZATION, DISINFECTION AND STERILIZATION

*Turn off the unit before any cleaning and disconnect the power cord from the mains.*

### CLEANING THE DEVICE AND THE TUBE'S OUTER SURFACE

*Only use a cloth dampened with antibacterial detergent (non-abrasive and free of solvents of any kind).*

## SANITIZATION, DISINFECTION AND STERILIZATION OF THE ACCESSORIES

We recommend a personal use of both the nebulizer and accessories to avoid any risk of infection or contagion. Open the nebulizer by turning the top (C1.3) counterclockwise, then remove the nozzle from the top (C1.3), and disassemble it as shown in the "Connection Diagram" section C1.

### SANITIZATION

Before and after each use, sanitize the nebulizer and accessories according to one of the methods described below.

**(Method A): The accessories that can be sanitized with the method A are** (C1.1-C1.2-C1.3-C1.4-C2-C3-C4) Sanitize the accessories in warm (approximately 40° C) drinking water with mild dishwashing liquid (non abrasive) or in the dishwasher with hot cycle

**(Method B): The accessories that can be sanitized with the method B are** (C1.1-C1.2-C1.3-C1.4-C2-C3-C4) Sanitize the accessories by immersion in a solution consisting of 60% water and 40% white vinegar. Finally rinse thoroughly with warm (approximately 40° C) drinking water.

**(Method C): The accessories that can be sanitized with the method C are** (C1.1-C1.2-C1.3-C1.4-C2) Sanitize the accessories by boiling in water for 20 minutes; preferably use of demineralized or distilled water to avoid scale deposits.

After sanitization shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).

### DISINFECTION

**The accessories that can be disinfected are** (C1.1-C1.2-C1.3-C1.4-C2-C3-C4)

The disinfection procedure described in this paragraph is to be performed before using the accessories and is effective only if all the prescriptions herein are met and all the components involved are previously sanitized.

The disinfectant to be used must be of electrolytic oxidizing chlorine type, (active ingredient: sodium hypochlorite), specific for disinfection and is available in all pharmacies.

#### Procedure:

- Fill a container, suitable to contain all the individual components to be disinfected, with a solution of drinking water and disinfectant, while respecting the proportions indicated on the packaging of the disinfectant itself.
- Completely immerse each individual component in the solution, taking care to avoid the formation of air bubbles in contact with the components. Leave the components immersed for the period of time indicated on the packaging of the disinfectant, and associated with the concentration chosen to prepare the solution.
- Remove the components now disinfected and rinse thoroughly with lukewarm drinking water.
- After disinfection shake the accessories vigorously and place them on a paper towel, or alternatively, dry them with a jet of hot air (e.g. using a hair dryer).
- Dispose of the solution as recommended by the disinfectant manufacturer.

### STERILIZATION

**The accessories that can be sterilized are** (C1.1-C1.2-C1.3-C1.4-C2)

The sterilization procedure described in this paragraph is effective only if all the prescriptions herein are met, all the components involved are previously sanitized and the procedure itself is validated according to ISO 17665-1 standard.

**Unit:** Steam sterilizer with fractionated vacuum and overpressure according to EN 13060 standard.

**Procedure:** Package each component to be treated in sterile barrier system or packaging complying with EN 11607 standard. Place the packed components in the steam sterilizer.


Run the sterilization cycle according to the instructions manual of the device, selecting a temperature of 134° C and a duration of 10 minutes

**Storage:** Store sterilized components as the instructions of the sterile barrier package or system chosen.

## AIR FILTERING

The unit is equipped with an air filter (A3) which should be replaced when dirty or when its colour changes. To replace the filter take it out by pulling it as shown in figure. The filter is designed to remain fixed in its housing.

**Use only Flaem original accessories.**

## SYMBOLS

	Class II device		Type BF applied part
	Important: check the operating instructions		Alternating current
	Switch off		Hazard: electrocution. Consequence: Death. Do not use device while taking a bath or shower
	Switch on		CE Medical Marking ref. Dir. 93/42 EEC and subsequent amendments
	Complies with: European standard EN 10993-1 "Biological Evaluation of Medical Devices" and European Directive 93/42/EEC "Medical Devices"		TÜV Certification ref. EN 60601-1 3° Ed. EN 60601-1-11
IP21	Enclosure protection rating: IP21. (Protected against solid bodies over 12 mm. Protected against access with a finger; Protected against vertically falling water drops.)		Serial number of device
	Keep dry		Manufacturer

## TECHNICAL FEATURES

### Aerosol therapy unit Mod. PO406EM F400

Voltage:	230V ~ 50Hz 130VA	115V ~ 60Hz	220V ~ 60Hz
In conformity with:			
Safety certifications:			

Power supply:	230V ~ 50Hz 130VA
Max pressure:	1.8 ± 0.3 bar
Compressor air output:	9 LPM approx
Sound level (at 1 m):	55 dB (A) approx
Operation:	continuous
Operating conditions:	Temperature: min 10°C; max 40°C RH Air humidity: min 10%; max 95% Atmospheric pressure: min 69KPa; max 106KPa
Storage conditions:	Temperature: min -25°C; max 70°C RH Air humidity: min 10%; max 95% Atmospheric pressure: min 69KPa; max 106KPa
Dimensions (W)x(D)x(H)	21x16x13 cm
Weight:	1.73 Kg

### Applied parts

Type BF applied parts are: patient accessories (C2,C3,C4)

### RF6 Basic<sup>2</sup> nebuliser

Medication minimum capacity:	2 ml
Medication maximum capacity:	8 ml
Operating pressure (with neb.):	0.60 bar

MMAD (µm) <sup>(2)</sup>		BREATHABLE FRACTION < 5 µm (FPF) <sup>(2)</sup>		DELIVERY ml/min <sup>(1)</sup> (approx)	
Standard mode	"High speed" mode	Standard mode	"High speed" mode	Standard mode	"High speed" mode
3.8	4.00	62%	63%	0.25	0.32

(1) data detected according to Flaem I29-P07.5 internal procedure

(2) In vitro characterization certified by Inamed Research GmbH & Co. KG on behalf of TÜV Rheinland Product Safety GmbH - Germany in compliance with the new European Standard for aerosol therapy units, Standard EN 13544-1, ANNEX CC. Further details are available on request.

**DISPOSAL OF DEVICE:** In conformity with Directive 2002/96/EC, the symbol shown on the device to be disposed of indicates that it is considered as waste and is therefore subject to “sorted waste collection”. The user must therefore take (or have taken) the above waste to a pre-sorted waste collection centre set up by the local authorities, or else give it back to the dealer when purchasing a new appliance of the same type.

Pre-sorted waste collection and the subsequent treatment, recovery and disposal operations favor the production of appliances made of recycled materials and limit the negative effects of any incorrect waste management on the environment and public health. The unlawful disposal of the product by the user could result in administrative fines as provided by art. 50 and subsequent amendments of Legislative Decree no. 22/1997.

#### **ELECTROMAGNETIC COMPATIBILITY**

This device has been designed to meet the current requirements for electromagnetic compatibility (EN 60 601-1-2:2007). Electromedical devices require special care during installation and use, with respect to EMC requirements. It is therefore required that they be installed and/or used according to the manufacturer's specification. Potential risk of electromagnetic interference with other devices, in particular with other devices for analysis and treatment Radio and mobile or portable RF communications devices (mobile phones or wireless connections) may interfere with the operation of electromedical devices. For additional information please visit [www.flaemnuova.it](http://www.flaemnuova.it). The Manufacturer reserves the right to make technical and functional changes to the product without notice.

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All rights reserved
**FLAEM NUOVA S.p.A.**