



AEROSOL SYSTEM COMPARISON CHART

	NASAL WASH	MR. BEETLE	MIZAR	BORA PLUS	MISTRAL PLUS	HELIOS	ALISEO PLUS
							
	Reciprocating piston compressor (designed for use without lubrication or maintenance)						Single chamber Ultrasonic Nebuliser
Power supply/ Power	3,0 V --- 4,5 w 2 Batteries AA (Alkaline/ Rechargeable)	230 V ~ 50 hz / 130 VA	12V DC 1,25 A	230 V ~ 50 hz / 140 VA	230 V ~ 50 hz / 130 VA	230 V ~ 50 hz / 130 VA	Power supply 100-240 V ~ 50/60 hz 15V DC Device 15 Vdc
Air filtering	-	✓	✓	✓	✓	✓	✓
Compressor Air output	-	~ 9 L/min	~ 9 L/min	~ 10 L/min	~ 9 L/min	~ 9 L/min	-
Max pressure	-	~ 1,5 Bar	~ 1,7 Bar	~ 2,2 Bar	~ 1,5 Bar	~ 1,5 Bar	
Sound level (a 1 m)	-	~ 55 Db (a)	~ 51 Db (a)	~ 57 Db (a)	~ 57 Db (a)	~ 54 Db (a)	< 50 Db (a)
Dimensions LxWxH (cm)	7,2x5,3x15,8	21x16x13	12x12x6	18x30x12	18x30x12	13,5x15,5x11	7x4x21
Weight	185 G (without batteries)	1,7 Kg	0,58 Kg	1,5 Kg	1,5 Kg	1,1 Kg	0,28 Kg
Transducer Frequency	-	-	-	-	-	-	2,4 Mhz
On/Off led	-						Green led light on No led light off
In Compliance With Dir. 93/42 CEE	-	CE 0051	CE 0051	CE 0051	CE 0051	CE 0051	CE 0051
In Compliance With Dir. 2006/95 CEE	CE	-	-	-	-	-	-
Approvals		 (Mod. P0406EM F400)		 (Mod. P0106EM F700)	 (Mod. P0106EM F400)	 (Mod. P0504EM F400)	

✓

YES

-

NO